

Specification

Disposable Arterial Blood Sampler				
				
Model	Size	Description	Needle Dimension	Packaging Quantity
isafeA-1	1mL	Pull-out type, include syringe protectors /syringe safety device/Luer Cap	25G*25mm	50pcs/pack 8packs/set
isafeA-3	3mL	Pull-out type, include syringe protectors /syringe safety device/Luer Cap	22G*25mm	50pcs/pack 8packs/set
isafeA-1N	1mL	Pull-out type, include syringe protectors /rubber safety tip/Luer Cap	23G*25mm	50pcs/pack 8packs/set
isafeA-3N	3mL	Pull-out type, include syringe protectors /rubber safety tip/Luer Cap	22G*25mm	50pcs/pack 8packs/set
isafeP-1	1mL	Preset type, include syringe protectors /syringe safety device/Luer Cap	23G*16mm	50pcs/pack 8packs/set
isafeP-3	3mL	Preset type, include syringe protectors /syringe safety device/Luer Cap	22G*25mm	50pcs/pack 8packs/set
isafeP-1N	1mL	Preset type, include syringe protectors /rubber safety tip/Luer Cap	23G*16mm	50pcs/pack 8packs/set
isafeP-3N	3mL	Preset type, include syringe protectors /rubber safety tip/Luer Cap	22G*25mm	50pcs/pack 8packs/set

Safer Choice Brings
the Best Result
Blood Sampler



Capillary Blood Sampler

Safer

- 100% plastic eliminates the risk of breakage for safer blood collection
- Replaces glass capillary tubes keep users safe

Accurate

- Dry-sprayed balanced heparin blood gas capillary tube avoids the dilution errors with liquid heparin
- High-concentrated lithium heparin minimizes the pre-analytical clots error

Variety

- Offers three sizes of capillary tubes: 85μL, 100μL, 115μL
- Mixing wire and magnet are available to meet different customers' need



Specification

Capillary Blood Sampler			
Model	Size	Description	Packaging Quantity
isafeC-85AN	85ul	Include capillary tube/cap	250pcs/pack 4 packs/set
isafeC-100AY	100ul	Include capillary tube/cap/mixing wire/magnet	250pcs/pack 4 packs/set
isafeC-100BN	100ul	Include capillary tube/cap	250pcs/pack 4 packs/set
isafeC-115BY	115ul	Include capillary tube/cap/mixing wire/magnet	250pcs/pack 4 packs/set

Disposable Arterial Blood Sampler

Reliable

- Ensure the accuracy of results by eliminating the interference of ions in the blood through precise ultrasonic spraying of calcium-balanced heparin-lithium.
- The conical part of the sampler can evenly dispense anticoagulant, reducing the risk of blood clotting.
- Luer lock provides good sealing to guarantee sample quality.

Security

- Equipped with safety needles to protect medical staff from accidental puncture.
- Screw connection to prevent needle from falling out.

Variety

- Two series of arterial blood sampler: isafeA (pull-out), isafeP (preset)
- Two sizes of arterial blood sampler tubes: 1mL, 3mL
- 1ml and 3ml are available with safety needle locker or rubber safety tip.

Suitable for all types of blood gas and chemistry analyzers, easy to operate in the clinic

