

# Rewrite the rules of physical therapy

Studies have shown that nearly 30 to 50 percent of women suffer from pelvic floor dysfunction, particularly older women and someone who has experienced childbirth. For the purpose of offering all-around solution for women care, the P series are launched into market to address the unmet need in pelvic muscle rehabilitation, which can provide precise assessment for pelvic muscle and effective physical treatment. It provides unmatched experience to both therapists and patients across a broad range of applications:

- Lower urinary tract symptom
- Disorder of defecation
- Pelvic organ prolapse
- Sexual dysfunction
- Chronic pelvic pain
- Diastasis Recti

## Available model

Model	Main function	Channels
PE4	Assessment	4 EMG/STIM Channels 1 Pressure Channel
P2	Assessment & Treatment	2 EMG/STIM Channels
P4	Assessment & Treatment	4 EMG/STIM Channels
P4 Pro	Assessment & Treatment	4 EMG/STIM Channels 1 Pressure Channel

## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



**Global Headquarters:**  
Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,  
Kengzi Subdistrict, Pingshan District, Shenzhen | 518122 P.R. China  
+86.755.26898326 | [www.edan.com](http://www.edan.com) | [info@edan.com.cn](mailto:info@edan.com.cn)

**U.S. and Canada inquiries:**  
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126  
+1.858.750.3066 | [www.edandiagnosics.com](http://www.edandiagnosics.com) | [edan-info@edandiagnosics.com](mailto:edan-info@edandiagnosics.com)

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

# A Moment of Intimacy

## P Series

Biofeedback and Stimulation System



ENG-OB&GYN  
P Series-V0.1-20191030

# P Series

## Biofeedback and Stimulation System

Biofeedback and stimulation system is a comprehensive solution for muscle performance assessment, re-education and rehabilitation. Aimed at offering more integrated clinical solutions for women health care, the P series is born to help physiotherapists to treat pelvic floor dysfunction for women of all age groups.




Electrode hook



Current control knob



Remote intensity control



## Muscular assessment

- **EMG & Pressure dual-mode evaluation**  
4 EMG/STIM channels and one pressure channel are provided for therapists.
- **Smart signal extraction algorithm**  
Intelligently extract pelvic muscular signal and exclude individual influence.
- **Scene simulation evaluation schemes**  
Lead patients to simulate different physiological activities for assessing muscular performance more comprehensively.

## Physical treatment

- **Intelligent therapeutic package**  
Generate recommended therapeutic package according to individual assessment result.
- **Flexible protocol adjustment**  
Support the therapist to adjust the protocols or parameters in real time.
- **Diverse treatment methods**  
Provide a variety of treatments to meet clinical needs, which contain electrostimulation, biofeedback training, relaxation training and dilation training.



Vaginal Stress Probe



SV1/SV2  
(Reusable/Disposable vaginal EMG /Stimulation sensor)



SA  
(Reusable rectal EMG/ Stimulation sensor)



Adhesive sensor  
(Disposable skin surface EMG/ Stimulation sensor)