

# F15 Air

## Comprehensive wireless monitoring solution

All fetal and maternal parameters can be monitored wirelessly, empowering mom and clinician with more choices, providing a more comfortable experience throughout the labor process.



# Superior Technology Incredible Performance

## F15

Fetal & Maternal Monitor



### About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



**Global Headquarters:**  
Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen  
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

**U.S. and Canada inquiries:**  
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126  
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com



© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

ENG-OB&GYN-F15-V1.1-20220725



# Evolutionary Generation of Fetal Monitoring

As an ultimate CTG machine, F15 series fetal and maternal monitor is the perfect combination of advanced technology and sophisticated algorithm for Obstetrics' new era. With state-of-art design, elegant appearance and strong mechanical structure, F15 Series is dedicated to bring the most efficient and innovative work flow to clinics.



15 mins Backup Memory



Transducer Plug and Play




Triplets Monitoring




<100g

## Sophisticated CTG Alarm System

The essential purpose of fetal monitoring is to recognize abnormal conditions in time. Integrated with international CTG interpretation standard, F15 will issue an alarm when dangerous FHR waveform is recognized.



Sinusoidal FHR



Low Variability FHR

## Innovative wireless application

Offering comprehensive wireless solutions to clinics, freeing midwives and mom from messy wires. Additionally, F15 Air supports non-stop monitoring when transferring mom-to-be, data can be transferred from one device to another seamlessly.



## FT20

Integrated maternal parameters module



- Support NIBP, SpO<sub>2</sub> and PR
- Wearable device provides mom with more choices with their ambulation and position



SpO<sub>2</sub>




NIBP



PR

## FECG-TA

As an advanced module for F15 Air, FECG-TA can acquire FHR and MHR from advanced surface electrodes system using fECG and mECG technology, and obtain EHG simultaneously. With all in one design, moms and caregivers are bestowed with more choices when monitoring.



Extraction

### Accurate


- Autonomous mECG channel for calibration
- Precise R-R interval identification from ECG
- Reliable for women of high BMI

### Comfort

- Cover all fetal positions
- Weight less than 150g
- Passive acquisition and non-invasive



Separation



Calculation

### All in one design

- fFECG, mECG, FHR, MHR, UA and MFM can be obtained with one module.