


Extraordinary maneuverability

HD colposcope are equipped with rolling stand or swing arm stand for smooth running and fine focusing. Rolling stand can be adjusted in vertical movements and swing arm is easier to move up and down for more stable and convenient positioning. Accurate and fast operation is just an easy case!

C6 HD video colposcope
One stop solution with integrated computer



C6A HD video colposcope
Flexible solution



About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Global Headquarters:
Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen 518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-OB & GYN-C6 HD&C6A HD
V1.1-20230206

C6 HD/C6A HD

Video Colposcope



C6A HD

C6 HD




HD Video Colposcope

This new HD digital colposcope serial is a perfect combination of outstanding HD camera, powerful software and convenient workstation to accommodate the rigid requirements of gynecologists in colposcopic screening. High resolution camera and unique LED light type ensure superior image quality to reproduce the examination area to the greatest extent. With professional software and analysis tool, everything becomes easy, standardized and well-regulated. Different models and accessories provides more available options to meet various demands in clinical application. This new serial is turning a new page in colposcopy.

High image quality to ensure flawless display at users' end

The outstanding full HD image quality contributes for the correct decision-making and provides the safety guarantee for the patient. Resolution is more than 900 TVL. HDMI video is also supportive for teaching and case discussion. Flexible magnification level and better vision scope provides extra clear view from different angels. High-precision technology provides superior image quality, improved clearness and focal range to facilitate the physicians' work.



Green light for the better visualization of tissue

Green light can enhance the contrast to highlight the abnormal blood vessels. The improved visualization can facilitate the physician to define the problem area.



Illumination

Our HD colposcope is equipped with 12 SMD LED in ring shape with adjustable intensity to provide ultra-bright LED illumination.



Patient management software

