

# Blood Sampler

Safer Choice Brings the Best Result



# Blood Sampler

## Safer

- 100% plastic eliminates the risk of breakage for safer blood collection
- Replaces glass capillary tubes keep users safe

## Accurate


- Dry-sprayed balanced heparin blood gas capillary tube avoids the dilution errors with liquid heparin
- High-concentrated lithium heparin minimizes the pre-analytical clots error

## Variety

- Offers three sizes of capillary tubes: 85µL, 100µL, 115µL
- Mixing wire and magnet are available to meet different customers' need

## Specification

Storage and Expiration	2°C - 40°C; 24 months
Sample volume	85µL, 100µL, 115µL
Packaging	250 tests/kit (mixing wire and magnet are optional)
Anti-coagulant	Lithium Heparin



## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



A world of potential

### Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen  
518122 P.R. China | +86.755.26898326 | [www.edan.com](http://www.edan.com) | [info@edan.com](mailto:info@edan.com)

### U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126  
+1.858.750.3066 | [www.edandiagnostics.com](http://www.edandiagnostics.com) | [edan-info@edandiagnostics.com](mailto:edan-info@edandiagnostics.com)

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-POCT-Blood Sampler  
V1.0-20230330