



iM3s
Vital Signs Monitor

Pocket spot checks at all times

iM3s

Vital Signs Monitor



About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



A world of potential

Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-iM3s(L)-V1.0-20221130



Small outside, smart inside

Being the smallest vital signs monitor among its kind, iM3s presents a fashionable handheld design, delivering extraordinary portability throughout spot checks. Though small in size, it integrates all necessary vital measurements together, including SpO₂, NIBP and optional TEMP.

A Dedicated Spot Monitor

With its pleasing size and large capacitive touch screen, iM3s could be operated by one hand, bringing comfortable and familiar user experience to daily applications. For each test, all measurements and custom parameters could be displayed on the same screen, while the test results being saved automatically in its built-in non-volatile memory or transfer to hospital information system.

 300g	 IP44	 1M	 ≥8 hours
Lightweight Design	Waterproof & Dustproof	Shock Resistance	Minimum Run Time
 Wi-Fi	 e-Link	 20,000 sets	
Network Connection	Wireless Connection	Internal Storage	Barcode Scanner



Tabletop Expansion

Plugged onto an optional tabletop dock, the iM3s easily fits into regular tabletop applications, enabling more possibilities, including AC charging, optional thermal recording, optional temporal temperature measurements, etc.



Powerful functions with flexible configuration

Main Unit

- Standard Configuration: SpO₂, NIBP, Capacitive touch screen, Li-ion battery
- Optional Configuration: Built-in barcode scanner, 2.4G/5G Hz Wi-Fi, e-Link

CS-04 Charger Stand

- AC power supply

CS-05 Charger Stand

- Standard Configuration: Built-in Li-ion battery, AC power supply, USB port, RJ45 port
- Optional Configuration: Thermal recorder, Exergen TEMP (TAT-5000S)