

# M3A Vital Signs Monitor



Outpatient



Transport



Ward



Portable Design



**EDAN**

A world of potential

# M3A Vital Signs Monitor

EDAN M3A Vital Signs Monitor provides essential measurement of patient's SpO<sub>2</sub>, NIBP and Temperature. It offers an effective solution and brings high performance vital signs monitoring at your finger tips.

Measured Parameters: NIBP + SpO<sub>2</sub> + Infrared Ear Thermometer

- 3.5" high resolution color TFT-LCD screen & LED display
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Real-time measurements display
- Trend table display
- Built-in rechargeable Lithium-ion battery for 18 hours of continuous work

**12,000**


Measurements  
Review

- Compatible with CMS-Lite for PC data management
- Nurse call
- Suitable for adult, pediatric and neonatal patients

## Seamless Connectivity



## Mounting Solutions



## Specification

Configuration	SpO <sub>2</sub>	NIBP	Infrared Ear Thermometer
SpO <sub>2</sub> only	✓	—	—
NIBP only	—	✓	—
SpO <sub>2</sub> +NIBP	✓	✓	—
NIBP+Infrared Ear Thermometer	—	✓	✓
SpO <sub>2</sub> +NIBP +Infrared Ear Thermometer	✓	✓	✓



### Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen  
518122 P.R. China | +86.755.26898326 | [www.edan.com](http://www.edan.com) | [info@edan.com](mailto:info@edan.com)

### U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126  
+1.858.750.3066 | [www.edandiagnostics.com](http://www.edandiagnostics.com) | [edan-info@edandiagnostics.com](mailto:edan-info@edandiagnostics.com)

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-M3A-IET-V1.9-20230516